

XD3 Hardware Specifications

Workstation

| Item | Minimum Specification | Recommended Specification |
|------------------|--|--|
| Processor | Intel i5 or equivalent | Intel Quad Core i7 or equivalent |
| Memory | 4GB RAM | 8 or 16GB RAM |
| Disk | 160GB. 10GB available for system. 7200 RPM | 500GB. 10GB available for system. 10k or 15k RPM |
| Graphics | 128MB graphics card, DVI | 512MB graphics card, dual DVI |
| Display | Single 1680x1050 flat panel | Dual 1680x1050 flat panel displays |
| Mouse | 2-button + wheel | 2-button + wheel |
| Keyboard | US or UK keyboard | US or UK keyboard |
| Operating system | Windows XP (32-bit) or Windows 7 (32-bit) | Windows 7 (64-bit) |

Laptop

| Item | Minimum Specification | Recommended Specification |
|------------------|---|--------------------------------------|
| Processor | Intel Core 2 Duo, Core i3 or equivalent | Intel Quad Core i7 or equivalent |
| Memory | 3GB RAM | 8GB RAM |
| Disk | 80GB. 10GB available for system. 7200 RPM | 250GB SSD. 10GB available for system |
| Graphics | 128MB in-built graphics | 256MB graphics card, DirectX9 |
| Display | 1440x900 | External display, 1680x1050 |
| Mouse | Laptop trackpad | External mouse, 2-button + wheel |
| Keyboard | Laptop keyboard | Docking station + external keyboard |
| Docking Station | | DVI external monitor support |
| Operating system | Windows XP (32-bit) or Windows 7 (32-bit) | Windows 7 (64-bit) |

Thin Client Server

Mirada will individually recommend server specifications to ensure suitability based on environment and number of concurrent users. This table is for guidance only.

| Item | Minimum Specification | Recommended Specification |
|------------------|--|--|
| Processor | Quad Core Xeon or equivalent, 2.50GHz | 2 x Quad Core Xeon or equivalent, 3.0GHz |
| Memory | 16GB RAM | 32GB RAM |
| Disk | 250GB. 20GB available for system. 7200 RPM | 750GB. 20GB available for system. 10K or 15K RPM |
| Graphics | 128MB integrated graphics | 128MB integrated graphics |
| Display | 1280x1024 flat panel for administrator | 1280x1024 flat panel for administrator |
| Mouse | 2-button | 2-button |
| Keyboard | US or UK keyboard | US or UK keyboard |
| Operating system | Windows Server 2008 R2, SP1 | Windows Server 2008 R2, SP1 |
| Network | 1Gb connection to client machines | 1Gb connection to client machines |