

XD3 Hardware Specifications

Workstation

Item	Minimum Specification	Recommended Specification
Processor	Intel i5 or equivalent	Intel Quad Core i7 or equivalent
Memory	4GB RAM	8 or 16GB RAM
Disk	160GB. 10GB available for system. 7200 RPM	500GB. 10GB available for system. 10k or 15k RPM
Graphics	128MB graphics card, DVI	512MB graphics card, dual DVI
Display	Single 1680x1050 flat panel	Dual 1680x1050 flat panel displays
Mouse	2-button + wheel	2-button + wheel
Keyboard	US or UK keyboard	US or UK keyboard
Operating system	Windows XP (32-bit) or Windows 7 (32-bit)	Windows 7 (64-bit)

Laptop

Item	Minimum Specification	Recommended Specification
Processor	Intel Core 2 Duo, Core i3 or equivalent	Intel Quad Core i7 or equivalent
Memory	3GB RAM	8GB RAM
Disk	80GB. 10GB available for system. 7200 RPM	250GB SSD. 10GB available for system
Graphics	128MB in-built graphics	256MB graphics card, DirectX9
Display	1440x900	External display, 1680x1050
Mouse	Laptop trackpad	External mouse, 2-button + wheel
Keyboard	Laptop keyboard	Docking station + external keyboard
Docking Station		DVI external monitor support
Operating system	Windows XP (32-bit) or Windows 7 (32-bit)	Windows 7 (64-bit)

Thin Client Server

Mirada will individually recommend server specifications to ensure suitability based on environment and number of concurrent users. This table is for guidance only.

Item	Minimum Specification	Recommended Specification
Processor	Quad Core Xeon or equivalent, 2.50GHz	2 x Quad Core Xeon or equivalent, 3.0GHz
Memory	16GB RAM	32GB RAM
Disk	250GB. 20GB available for system. 7200 RPM	750GB. 20GB available for system. 10K or 15K RPM
Graphics	128MB integrated graphics	128MB integrated graphics
Display	1280x1024 flat panel for administrator	1280x1024 flat panel for administrator
Mouse	2-button	2-button
Keyboard	US or UK keyboard	US or UK keyboard
Operating system	Windows Server 2008 R2, SP1	Windows Server 2008 R2, SP1
Network	1Gb connection to client machines	1Gb connection to client machines