

Casemeeting Hardware Specifications

Casemeeting Server

Mirada will individually recommend server specifications to ensure suitability based on environment and number of concurrent meetings. This table is for guidance only.

Item	Minimum Specification	Recommended Specification
Processor	Quad-core Xeon @ 2.5GHz or equivalent	2 x Quad-core Xeon @ 3.0GHz or equivalent
Memory	8GB RAM	24GB RAM
Disk	500GB	1TB to 2TB (RAID-1 or RAID-5)
Operating system	Windows Server 2008 R2 SP1	Windows Server 2008 R2 SP1

Casemeeting uses RDP to view the presentation. A fast network is recommended, with an end-to-end connection speed of 1Gb and low latency.