

Cedars-Sinai Cardiac Suite Hardware Specifications

The following minimum software and hardware requirements should be satisfied before installing CSMC Cardiac Suite.

Function	Specification
Operating System	Windows 7 (32 & 64 bit): Home Premium, Professional, enterprise and Ultimate. Windows Server 2008 R2 (32 & 64 bit): Standard Edition and Enterprise Edition. Windows Server 2008 (32 & 64 bit): Standard Edition and Enterprise Edition. Windows Vista (32 & 64 bit): Home Premium, Business, Enterprise and Ultimate. Windows Server 2003: Standard Edition and Enterprise Edition. Windows XP: Home, Professional 64-bit. Mac OS X: 10.6 (Snow Leopard) and 10.7 (Lion).
Ram (Random Access Memory)	1 GB (2 GB for Fusion/CT)
Available Disk Space	150MB for installation, additional space required for image data storage
Screen Resolution	1024x768 with 16-bit color. Wide-screen displays that meet the minimum requirements are supported.
Network Port	Ethernet network adapter (only required for workstation networking scenarios)
Miscellaneous	Mouse (or other pointer devices such as trackpad, trackball, etc) Keyboard