

RTx™ 1.8 Hardware Specifications

WorkStation

Item	Recommended Specifications
Processor	Intel Quad Core i7 or equivalent
Memory	16GB RAM
Disk	500 GB. 10GB available for system 10k or 15k RPM
Graphics	512MB graphics card, dual DVI
Display	Dual 1680x1050 flat panel displays
Mouse	2-button + wheel
Keyboard	US or UK keyboard
Operating System	Windows 7 (64-bit) or Windows 8.1 (64-bit) or Windows 10 (64-bit) Both physical workstations and virtual machines are supported

Thin Client Server

Mirada will individually recommend server specifications to ensure suitability based on environment and number of concurrent users. This table is for guidance only. Client machines should be reasonably modern but do not require a high specification and must run Windows 7 (SP1 recommended), Windows 8.1, or Windows 10.

Item	Recommended Specifications
Processor	Three cores/2.5GHz per concurrent user
Memory	16GB RAM per concurrent user
Disk	750GB. 20GB available for system 10k or 15k RPM
Graphics	128MB dedicated graphics
Operating System	Windows Server 2008 R2, SP1 Windows Server 2012 R2 Both physical servers and virtual machines are supported
Network	1Gb connection to client machines
Network Card	Dual teamed cards for four or more concurrent users



mirada-medical.com

Worldwide Headquarters
Mirada Medical Ltd.
Oxford Centre for Innovation
New Road, Oxford, OX1 1BY
United Kingdom

info@mirada-medical.com
+44 (0)1865 261410

USA Office
Mirada Medical USA, Inc.
999 18th Street, 2230S
Denver, CO 80202
United States

info@mirada-medical.com
+1 877.872.2617