

Invia Corridor 4DM Hardware Specifications

Item	Specifications
Operating System	XP, Vista, Windows 7, Windows 8- applicable to both 32 and 64 bit
Processor Speed	Dual Core 2.0 GHz or greater
Bus Interface	PCI Express x16 or x1 depending on the model
PCI Express	2.0 Support
OpenGL	Video Card must support OpenGL 3.0 greater (4.1 preferred)
Internal System Memory	2 GB minimum, 3GB recommended
Video Memory	256 MB for SPECT; 512 MB or more for PET and CT
Monitor Display	24-bit color; 1280 X 1024 recommended; 20-inch or greater
Video Card Max Digital Display	1920 X 1200 (DVI); 2560 X 1600 (Display Port)
Video Card Max Analog Display	2048 X 1536
Disk Space	1 GB required for installation; 2-4 GB for image files



mirada-medical.com

Worldwide Headquarters
Mirada Medical Ltd.
Oxford Centre for Innovation
New Road, Oxford, OX1 1BY
United Kingdom

info@mirada-medical.com
+44 (0)1865 261410

USA Office

Mirada Medical USA, Inc.
999 18th Street, 2230S
Denver, CO 80202
United States

info@mirada-medical.com
877.872.2617