

Workflow Box 2.0 Hardware Specifications

Workflow Box Server (without DLCExpert)

Item	Recommended Specifications
Processor	2 x Quad-core Xeon or equivalent
Memory	32GB RAM
Disk	250GB SSD 20GB available for system
Operating System	Windows Server 2012 R2 Windows Server 2016 Both physical workstations and virtual machines are supported

Workflow Box Server (with DLCExpert)

In order to optimise machine learning functions, a non-integrated graphics card is highly recommended.

Item	Recommended Specifications
Processor	2 x Quad-core Xeon or equivalent
Memory	32GB RAM
Disk	250GB SSD 20GB available for system
Graphics Card	GeForce GTX 1080, equivalent or greater Graphics memory of 6GB or greater CUDA compute Capability: 6.1
Operating System	Windows Server 2016

Note: The following graphics cards have been tested by Mirada Medical and confirmed to work with DLCExpert:

- NVIDIA GeForce GTX 1060 6GB
- NVIDIA GeForce GTX 1080 8GB
- NVIDIA GeForce GTX 1080Ti



mirada-medical.com

Worldwide Headquarters
Mirada Medical Ltd.
Oxford Centre for Innovation
New Road, Oxford, OX1 1BY
United Kingdom

info@mirada-medical.com
+44 (0)1865 261410

USA Office

Mirada Medical USA, Inc.
999 18th Street, 2230S
Denver, CO 80202
United States

info@mirada-medical.com
+1 877.872.2617